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Fax +49 74 31 9 89 75 15

<http://www.ergosana.de/>

Couch/Semi-Couch Safety Ergometer Sana couch 150 / 250 L and  
Stress Echo Safety Couch Sana cardio 150 / 250 SE

---

## User Guide


**CE Conformity Declaration**

ergosana GmbH herein declares that the products of the ergometer system Sana couch 150 / 250 L and Sana cardio 150 / 250 SE have been designed and manufactured in accordance with the relevant requirements of the EC directive 93/42/EEC.

This declaration loses its validity if the above devices are modified without ergosana's consent.

Manufactured by:  
ergosana GmbH  
Truchtelfinger Str. 17  
D-72475 Bitz

CE0124

.....  
Harald Neukirchner  
Quality manager

## Table of Contents

<b>1 General</b>	<b>4</b>
1.1 Intended use	4
1.2 Contra-indications	4
1.3 Features	4
1.4 Instruction	4
1.5 Maintenance	4
<b>2 Product description</b>	<b>5</b>
2.1 Unit Components	5
2.2 Accessories	5
2.3 Potential Equalisation	5
2.4 Technical Data	6
2.5 Signs and Symbols	6
<b>3 Installation</b>	<b>7</b>
3.1 Location	7
3.2 Assembly Instructions	7
3.2.1 Unpack and Assemble the Running Gear	7
3.2.2 Mount Blood Pressure Recorder and Armrest for 250 SE	7
3.2.3 Set Saddle and Headrest	7
3.2.4 Connect Device	7
<b>4 Unit Components</b>	<b>8</b>
4.1 Control Panel	8
4.2 Rotational Speed Display on the Control Console	8
4.2.1 Panel Interfaces for 250 SE	8
4.2.2 Keys and Display	9
4.2.3 Setting the Language	10
4.3 Blood Pressure Cuff for 250 SE	10
4.4 Remote Control	11
5.2 Safety Precautions when Operating with Other Devices	13
5.3 Maintenance Safety Precautions	13
5.4 Interference	13
<b>6. Initial Operation</b>	<b>14</b>
6.1 Blood Pressure Measuring Unit (250 SE)	14
6.2 Applying the Cuff (250 SE)	15
6.3 Positioning of the Patient on the Couch	15
<b>7. Ergometry</b>	<b>16</b>
7.1 Defining the Automatic Load Programs	16
7.2 Recommended Settings	18
7.3 Remote Operation	19
7.3.1 Explanation	19
7.3.2 Prerequisites	19
7.4 Training Program (Option)	19
7.4.1 Training with Constant Heart Rate (Pulse-Steady-State)	20
7.4.2 Configuring a Training Program on the Ergometer	20
<b>8 Maintenance and Fault Clearing</b>	<b>22</b>
8.1 Measurement Check	22
8.2 Cleaning the Device	22
8.3 Cleaning the Blood Pressure Cuff (250 SE)	22
8.4 Repairing Electromagnetic Disorders	23
8.5 Checking and Setting the Supply Voltage	24
8.6 Changing a Mains Fuse	24
8.7 Disposal Information	24
<b>9 Appendix</b>	<b>25</b>
9.1 Technical Customer Service and Sales Locations	25

## 1 General

The couch/semi-couch safety ergometer *Sana couch 150 L / 250 L* and the stress echo couch *Sana cardio 150 SE / 250 SE* are state-of-the-art high performance ergometers. The models designated Sana couch 250 L and Sana cardio 250 SE are additionally equipped with a blood pressure measuring module. This module is situated in the ergometer's control console. The device meets the highest quality standards in precise physical exertion tests for conducting measurements in the field of cardiovascular and pulmonary function diagnostics. High safety standards are ensured thanks to the patient's reclined position.

### 1.1 Intended use

Ergometer Sana bikes are bicycle ergometers that are intended for defined exercise ergometry during a patient's examination and therapy. These products are used in practices, clinics, therapy and rehabilitation centres. The ergometers are operated by physicians and medical personnel.

### 1.2 Contra-indications

**In the event of the following contra-indications, NO exercise test must be performed:** in the case of acute cardiac infarction or unstable angina pectoris, serious hypertonia at rest, carditis, insufficiency of the heart, serious valvular heart defect, serious cardiac arrhythmia at rest, aortic aneurysm or other manifest cardiovascular diseases.

### 1.3 Features

**The following characteristics make the unit exceptional:**

- Attractive design
- Easy access supported by a footboard and safe handlebar
- Stable steel construction, compact drive unit
- Comfortable, stable couch enabling the patient to lie still
- Safe position of the patient in the case that resuscitation should be necessary
- Easy adjustment of couch, saddle and headrest using a handy remote control
- New, high-performance control electronics
- Graphic display featuring visual representation of ergometry data
- Easy operation via menu guidance
- Remote operation – personalised programs – training programs
- Disturbance-free blood pressure measurement
- Performance range from 1 to 999 watts
- Guaranteed accuracy (error factor < 3% in the independent rpm range)
- Very quiet drive unit
- Pleasant pedalling sensation due to large gyrating mass
- Galvanically isolated RS-232 interface, secure data transfer

### 1.4 Instruction

Before the initial operation, carefully read through this user guide, paying special attention to the warnings and safety instructions.


### 1.5 Maintenance

This is a low-maintenance device. You will find detailed maintenance instructions in section 8.

## 2 Product description

### 2.1 Unit Components

1. Paper roll
2. Blood pressure cuff for 250 SE
3. Variable armrest for 250 SE
4. Remote control
5. Couch
6. Footboard
7. Running gear with lockable heavy rollers
8. Blood pressure unit connection for 250 L/SE
9. Control panel with LCD display and foil keyboard with control elements and LED rotational speed display
10. Handlebar as mounting aid
11. Handlebar
12. Adjustable saddle
13. Adjustable headrest with paper roll carrier

ENGLISH

### 2.2 Accessories

#### Every device comes with:

- Power cable with European plug
- Blood pressure cuff for 250 SE
- Armrest for blood pressure measurement for 250 SE
- Handlebar
- Saddle
- Headrest with paper roll carrier
- Paper roll
- User guide
- Inspection report

### 2.3 Potential Equalisation

A standard potential equalisation stud is located on the rear side, next to the power connection unit. It is marked with a green/yellow information sign. Using an earthing cable, the ergometer can be connected with the potential equalisation of the examining room, which serves as a common earthing point for all other mains-operated devices in the room to ensure that all devices have the same earthing potential.

## 2.4 Technical Data

Bicycle ergometer with blood pressure measurement in accordance with DIN 13405/DIN VDE 0750-238.

<b>Braking principle</b>	Computer-controlled brakes with permanent measurement of torque. Braking performance is independent of revolutions per minute. (see load range precision)
<b>Power range</b>	1 till 999 Watts
<b>Load range</b>	independent rpm range 20 till 999 Watts
<b>Range of revolutions</b>	30 to 130/min
<b>Load precision</b>	3%, not less than 3 watts (in the independent rpm range)
<b>Load parameters</b>	1. In keeping with set internal load program 2. Parameters from external master unit over interface, in 1 watt steps. 3. Manual in 5-watt and 25-watt steps
<b>Load software</b>	5 freely programmable ergometry programs
<b>Time intervals</b>	1 min to 99 min
<b>Display</b>	Graphic LCD with 320 x 240 pixels, CCFT background lighting
<b>Blood pressure measurement for 250 SE</b>	Indirect with special modified R-R measurement system, computer evaluation and distortion-free suppression of interferences during ergometry. Automatic pressure release by 3 mmHg/pulse; quick pressure release at the average of high amplitudes. Measuring range 40...300 mmHg.
<b>Pulse measurement</b>	Priority principle; 1st ECG, 2nd RR; measurement range 35 to 240 bpm; accuracy $\leq 2\%$ .
<b>Height of saddle and headrest</b>	Infinitely adjustable for heights between 120 cm and 210 cm; weight max. 160 kg
<b>Long-term accuracy</b>	Torque equalisation at any time with weight
<b>Power supply</b>	230 VAC 50–60 Hz, 115 VAC 50–60 Hz The unit is suitable for use in networks according to CISPR 11, group 1, class B.
<b>Electric inputs/outputs</b>	RS-232 (galvanically isolated), analogue remote boot
<b>Couch dimensions (horizontal)</b>	55 x 160 cm
<b>Weight</b>	102 kg

## 2.5 Signs and Symbols

In this section, the signs and symbols used in connection with this device are explained:

Mains operated, alternating current

Potential equalisation connection (earth)

BF classified component

Warning! Follow the instructions in the documentation.

93/42/EEC for medical products 0124  
DEKRA

Protection class of the casing: IPX0

## 3 Installation

### 3.1 Location

Install the device in a suitable position (refer to safety instructions in section 5).

The unit should not be stored or operated in wet, moist or dusty surroundings. Nor should the unit be exposed to direct sunlight or other sources of warmth.

The unit should not come into contact with acidic vapours or fluids.

The unit should not be placed near X-ray units, large transformers or electrical motors. There must be a distance of at least one meter between the unit and the mains network.

### 3.2 Assembly Instructions

#### 3.2.1 Unpack and Assemble the Running Gear

After unpacking the couch ergometer and unscrewing it from the transport container, turn the frame on the container so that the enclosed heavy rollers can be mounted on the running gear.

The device is delivered with the control console facing the front so that the display can be seen by the person operating the machine.

#### 3.2.2 Mount Blood Pressure Recorder and Armrest for 250 SE

Mount the armrest for blood pressure measuring onto the left and the handlebar onto the right guide rail. Connect the blood pressure cuff by screwing the hose connection onto the joining nipple in the control console and inserting the plug into the socket next to it. The cuff can be placed on the left or the right arm from this position. For measurements on the right arm, the armrest has to be mounted on the right guide rail.

#### 3.2.3 Set Saddle and Headrest

There are two bore holes in the saddle's and the headrest's guide tubes. With these, the saddle and the headrest can very easily be adjusted to the sizes of adults (up to 210 cm) and children (down to 120 cm). For adults' size, the saddle and the headrest are fixed through the lower of each pair of bore holes. In addition, the headrest can be adjusted manually to find a comfortable position for very small or very tall patients.

Couch, saddle and headrest can be varied via the remote control for different sizes, according to the imprinted symbols.

#### **Caution**

When adjusting the saddle to a higher position, it is recommended to relieve the hoist motor by shifting the body weight onto the pedals.

#### 3.2.4 Connect Device

Establish potential equalisation (see section 2.3) and plug the supplied power cable into an earthed socket. Switch the device on using the main switch on the rear side.

The couch ergometer is now ready for operation.

## 4 Unit Components

### 4.1 Control Panel

The control panel is mounted on the upper side of the control panel carrier with two plug-in tabs. During normal operation, the display should be facing the examiner. It is possible to turn the console 180 degrees for special applications such as patient training, etc., so that the patient can reach the control elements and see the display.

All of the control electronics for ergometer operation and blood pressure measurement are located in the control panel.

A backlit LCD featuring the complete range of information is located on the front side under a foil keyboard with a see-through window. The control elements for the adjustment and operation of the ergometer are located on the foil keyboard.


The connections for the blood pressure cuff are located on the underside of the 250 L/SE. In couch ergometers, they contain an extension line. The cuff socket is located on the couch's upper edge.

An LCD on which the patient can see the number of pedal rotations per minute is located on the panel's upper side.

### 4.2 Rotational Speed Display on the Control Console

$n$  = crank rotations per min


$U \text{ min}^{-1}$

#### 4.2.1 Panel Interfaces for 250 SE

The extension line for the blood pressure cuff is connected at the bottom of the control panel.

1. Cuff connection
2. Microphone connection


#### 4.2.2 Keys and Display

The LCD with graphics capability, 320 x 240 pixels and a surface of 100 mm x 75 mm is located on the front of the control panel. It is covered by a foil keyboard with a see-through window.

The display shows all configuration and operating procedures. During the exercise testing, all current measurement data is shown in the display in alphanumeric and graphic forms. It is thus especially easy to monitor exercise testing.

The cursor can be moved through the selection menu in the display with the up and down arrow keys.

The right and left buttons on the control panel have been arranged so that the function executed by the corresponding key is shown directly above the key on the lower side of the display.

- 1 = Blood pressure measurement key
- 2 = Stop/quick pressure release key
- 3 = Confirmation key
- 4 = "Up" and "down" cursor keys
- 5 = LCD display

### 4.2.3 Setting the Language

The preset menu language is "deutsch".

To change the language, highlight the "Einstellungen" menu option using the "ab" arrow, and confirm with "OK".

In the "Einstellungen" menu, select the "Sprache" option and confirm with "OK". Then select the desired language using the "auf" or "ab" arrow, and confirm with "OK".

The menu is now displayed in the selected language.

### 4.3 Blood Pressure Cuff for 250 SE

The standard blood pressure cuff (order no. 24-10-301) has Velcro fixing. It can be used for arms up to 40 cm in diameter. A larger cuff (order no. 24-10-321) is available for larger arm diameters.

A microphone is installed in a microphone pocket on the inside of the cuff. It serves the transmission of the blood pressure sound.

The connection cable with an air and a microphone connection is 110 cm long. This is a sufficient length. It ensures that the cable does not dash against the ergometer during pedalling. This is to prevent unnecessary artefacts that might result in inaccurate blood pressure measurement. Longer cables (200 cm) are available but it is imperative to make sure that artefacts are prevented.

#### **Cleaning:**

The cuff should only be washed with soapy water and immediately dried off again. The surface of the microphone pocket is waterproof. Make sure that no moisture penetrates the opening of the microphone pocket. If this occurs repeatedly, it may result in damage to the microphone.

#### 4.4 Remote Control


For the quick and easy variation of the couch, the saddle and headrest are operated using a push-button remote control. With this, the correct position for the patient can be found very easily.

The top two buttons (1) control the couch's slope.  
 Left button: Couch is elevated  
 Right button: Couch is lowered

The two middle buttons (2) control the headrest.  
 Left button: Headrest is moved downwards  
 Right button: Headrest is moved upwards

The two lowest buttons (3) control the saddle.  
 Left button: Saddle is moved downwards  
 Right button: Saddle is moved upwards

**Note – with Stress Echo Couch only**  
 These two buttons control the couch's tilt to the left.

ENGLISH

#### 4.5 Additional Instructions for the Stress Echo Couch

In contrast to an ordinary couch ergometer, the couch of the stress echo ergometer can be tilted infinitely by up to 45 degrees to the left. By this, the heart is brought into a suitable position for ultra-sound examinations during exercise. In this position, clearer and sharper ultra-sound images of the heart are achieved. The tilting is controlled via the remote control.

The following supports ensure the patient's safety in any tilt position and give him a secure feeling.

##### Lap Belt:

The belt is lead around the patient's lap and fixed using the Velcro strap in a position to keep the patient in the middle of the couch surface. On applying the belt, make sure that it does not hinder the patient's pedalling. When the belt is not used, roll it up and fix it to the ergometer's right side using a Velcro strap.

##### Handlebar with Shoulder Support:

This handlebar is positioned to keep the patient in the middle of the tilted couch. The handle is extended to a distance that keeps the patient's arm away from the examination area. The shoulder's angle of rotation must be considered.

##### Right Armrest

This armrest is set to support the patient's right arm in a comfortable position. For a stress echo examination, the blood pressure cuff must be applied on the right arm. Please make sure that the microphone is positioned at the inner arm on the brachial artery and the cuff is lying flat on the arm.

##### Blood Pressure Measurement during a Stress Echo Examination

For normal ergometry examinations, we recommend to fix the armrest to the left and the handle to the right. The blood pressure can be measured on the left arm as usual.

During a stress echo examination, the blood pressure must be measured on the right arm. For this, the cuff is applied so that the microphone is positioned on the right inner arm approx. 3 cm above the elbow. **The cuff must lie flat with the tube exit facing downward.** This is slightly more difficult on the right arm than on the left arm. This is why a little bit more care must be applied with this method. Further more, you must make sure that the cuff's air tube cannot dash anywhere.

##### Pedals and Pedal Straps

The feet must be positioned with the balls on the pedals. They are then secured using straps. Pedals in which the whole foot is entered and fixed are available as an option.

##### Couch's Installation and Stability

Position the couch in a way giving you access to the patient from any side. The couch's left back foot is extended and contains a base adjustor. Block all 3 rolls with the brakes. Adjust the base adjustor to prevent any play. This ensures the ergometer's absolute stability in all tilt positions up to a patient weight of 160 kg. If the couch should still shake due to floor unevenness, all rolls can be replaced by base adjustors.

## 5. Safety Notes

### 5.1 Precautions during Operation

Before using the unit, make sure that the medical product consultant has conducted an introduction in regard to function and safety precautions.

The unit should not be used if there are any doubts as to its being isolated from earth or the suitability of the power cable.

The supplied power cable meets valid regulations for medical applications.

The unit is not intended for use in wet rooms, outdoors or in areas where there is danger of explosion.

Before beginning the operation, the running gear's heavy rollers must be blocked by pressing down the brakes. Make sure that the unit is set up on flat ground to ensure absolute stability.

The ergometer must be located in a way to ensure that there is still enough room for operation when the couch is in a horizontal position and the headrest is fully extended.

The holding straps on the pedals must fit perfectly across the upper side of the shoe and be fastened with a Velcro strap.

### 5.2 Safety Precautions when Operating with Other Devices

When several devices are coupled, there is a risk that the leakage currents may add up.

The RS-232 interface, which can be used for communication with other devices, is isolated to ensure the patient's safety.

External devices may only be connected with the interface cables supplied by ergosana.

Portable communication devices, HF radios and devices labelled with the  symbol (non-ionic electromagnetic radiation) can affect the operation of this device (see section 8.4).

### 5.3 Maintenance Safety Precautions

The device must be turned off and the power plug disconnected before cleaning with liquid cleaning agents.

Use only standard cleaning agents for plastic and leatherette surfaces.

The unit may only be opened, repaired and serviced by authorised and trained personnel.

### 5.4 Interference

The unit meets EMC regulations for medical products to ensure protection against emission and radiation. Special caution should be taken when using this unit in combination with high-frequency devices.

## 6. Initial Operation

### 6.1 Blood Pressure Measuring Unit (250 SE)

In order to conduct exercise testing correctly, it is of utmost importance to measure physical performance data and data from the ECG measurement as well as simultaneously measuring and recording blood pressure data to determine the reaction of the circulatory system to increased exertion.

For this purpose, ergosana has developed an extremely accurate blood pressure measuring system that is not susceptible to interference. It has been integrated into this ergometer and uses a so-called indirect method of blood pressure measurement. The Korotkoff sound, which is created by the air being forced out of the cuff as blood flows through the area of compression, is recorded along with several other important parameters of critical importance to attain accurate measurement. These measurements are evaluated in milliseconds by an internal digital evaluation system and shown in the ergometer's display as systole and diastole. The pulse rate is also determined during measurement and likewise shown in the display. At the same time as they are shown in the display, the measurements can also be transferred to a peripheral device such as an ECG or pulmonary function unit for evaluation and recording over an RS-232 interface.

The blood pressure cuff is the measurement sensor for blood pressure. Despite the perfectly functioning measuring system, it remains critically important that the cuff is placed on the arm correctly and carefully. According to international agreement, the blood pressure should be measured on the left arm, which is near the heart, as the flow impedance level is lowest there. An exception to this rule is formed by approx. 1 to 2 per cent of test persons on whom the Korotkoff sound cannot be measured due to vascular phenomena. The cuff is placed on the right arm of such patients.

An armrest is attached to the ergometer for the arm on which the blood pressure is measured. The armrest provides the arm's stable position to prevent distortion of the result caused by artefacts. It must be set up to support the arm in a slightly angled position. It can be mounted on the right or the left side.

Please note that the cuff's air tube must be fixed in a way that prevents it from encumbering the ergometer. This is to prevent unnecessary artefacts that might affect the measurement's accuracy.

## 6.2 Applying the Cuff (250 SE)

The microphone is positioned so that it lies on the brachial artery, the largest arm artery. The location of the microphone in the cuff is marked with a red fabric tag.

The ideal location for the microphone is approx. 2 centimetres above the elbow joint on the inside of the arm, below the biceps. The cuff must be put on so it is tight and cannot shift out of position during the movement created during the stress test.

The cuff is inflated rapidly at the start of the measurement. The blood pressure and pulse rate are already roughly measured during pumping and the inflation pressure is determined.

After the systolic pressure value has been attained, the air is released from the cuff at a rate of 3 mmHg per heartbeat.

This procedure guarantees approximately equal measuring times despite the rising pulse rate during exertion.

The blood pressure measurement should not exceed a maximum total length of 45 seconds. One minute is indicated as the shortest measurement interval. A measurement interval of 2 or 3 minutes is preferable in most cases.

The blood pressure values are shown in the ergometer's graphic display together with the load and pulse rate graphs.

## 6.3 Positioning of the Patient on the Couch

The advantage of semi-couch ergometry is the patient's ideal positioning on the ergometer. In a 45 degrees position, the performance corresponds to the performance on a bicycle ergometer. Nevertheless, the patient is lying in a very relaxed position. This prevents artefacts during blood pressure and ECG measurement, enabling a better evaluation of the results. For elderly and sick patients, the reclined position is very safe and comfortable. The 45 degrees position should only be left if resuscitation is necessary or for heart catheter examination.


The saddle serves for setting the correct distance between the body and the pedals. The headrest's purpose is to provide a comfortable position for the head. The paper on the paper roll acts as surface cover and sweat absorber during the stress test.

## 7. Ergometry

This section describes ergometry with the internal ergometer software.

### 7.1 Defining the Automatic Load Programs


The following text field appears in the ergometer after it is switched on:

The menu item "Start ergometry" is highlighted with a black bar. This means the item is activated. It is possible to select an exercise program directly from here and begin ergometry.


Beforehand, however, the five varying exercise programs must be defined according to the examiner's needs and wishes. On delivery of the unit, each program is set with normal values which cannot harm the patient in case they are accidentally activated.

Activate "Define program" via the arrow keys  .


Confirm with the right "ok" button.

The program selection menu appears.

The exercise programs 1–5 are set according to the parameters of standardised ergometry or according to individual parameters in such a manner that the correct programs for the different patient groups can be called up at the touch of a button.

Confirm program 1 with the right button “ok”, and the window containing the settings menu for ergometry program no. 1 appears.

Check the existing settings for:

- Type of load (Stage or Ramp)
- Initial load
- Increase in load
- Stage time
- Blood pressure measurement interval (for 250 L/SE)
- Recovery load
- RR interval 2 min
- RR interval on

Push the “change” button if the settings need to be altered.

The individual parameters then appear one after the other. Select the desired values with the arrow buttons and confirm by pressing the “ok” button. The next parameter then appears. The settings menu can be exited at any time by pressing the “cancel” button.

## 7.2 Recommended Settings

An unlimited number of values can be saved in the programs 1 to 5, as shown in the examples below:

Program number	Initial load [watts]	Increase in load [watts]	Stage time [min]	BP interval [min]	Recovery load [watts]
1	30	10	1	2	20
2	25	25	2	2	25
3	50	25	2	2	25
4	50	50	3	3	50
5	75	50	3	3	50

The program always returns to the initial configuration after the individual programs have been set. "Start ergometry" is highlighted black on the LCD display. An ergometry exercise can be started immediately by pressing the "ok" button.

## 7.3 Remote Operation

### 7.3.1 Explanation

Remote operation means that the ergometer is externally controlled via the digital RS-232 interface or USB connector, i.e. all commands for the load and blood pressure measurement intervals are transmitted from a separate "master unit".

This type of operation is used primarily when several ergometers are controlled from a central workstation, the ECG device has its own ergometry-control program and the ergometer and ECG unit – and possibly other devices such as spirometry devices – can be combined to form an ergometry or spirometry station.

In combination with SCHILLER ECG units, remote operation is the only approved operation method.

### 7.3.2 Prerequisites

When the remote operation mode is selected, the unit used (an ECG device or a PC) must be connected with the ergometer via an interface cable. For our ergometers, this is a type RS-232 interface, which is isolated to ensure the patient's safety.

The appropriate baud rate must be selected under menu option "Setting > Interface". The appropriate transmission protocol must then be selected in the same menu under "Instruction Set". The so-called "ergoline operating mode" is stored under the setting P 10. Our own transmission protocol is located under "ergosana".

Information on interface and instruction set can be obtained from the operating data of the master unit. When the settings have been correctly defined, the ergometer automatically switches to "remote operation" when the first command comes through over the interface. The ergometry screen is opened in the display, showing the current load, blood pressure and pulse rate data in both alphanumeric and graphic forms. The execution of the internal programs is disabled in this operating mode.

Remote operating mode is ended by either pressing the command button "End" or by switching the device off.

## 7.4 Training Program (Option)

The use of the training program requires a pulse signal receiver (Polar system), which is integrated in the control panel. The receiver can be ordered together with the device or added later on.

The patient wears a transmitter belt, which is positioned on the skin below the chest. The effective radius of the signals transmitted by the belt and received in the control panel is approx. 70 cm. Please note that if the skin is dry, contact problems may occur between the belt and the skin at the beginning of the training. Therefore, if the pulse transmission is unsteady or disturbed, moisten the contact surfaces of the transmitter belt with contact spray or water.

#### 7.4.1 Training with Constant Heart Rate (Pulse-Steady-State)

A bicycle ergometer training with constant heart rate (pulse-steady-state method) in the individual training range is a highly efficient and risk-free training method for the cardiopulmonary system. Ask your physician which is the suitable training heart rate for you. He or she will determine it for you using an exercise test.

#### 7.4.2 Configuring a Training Program on the Ergometer

To define a training program, select the menu option **Define program**. Then select the menu option **Training** using the arrow keys. The below settings panel is displayed. To set or change the parameters in order, press **modify**.

The first setting defines the *Initial load*, e.g. 50 watts.

The setting *Duration A1* controls the duration of warming-up phase 1, e.g. 1min. The setting *Increase in load* controls by how many watts per minute the load is increased. This is warming-up phase 2.

In special cases, this phase can be limited in time using the next setting *Duration A2*. By default, the warming-up phase 2 is finished when the target HR is reached.

When the target HR is reached, the load applied is automatically reduced by 10%. This measure is to prevent a too great rate increase during the training phase which now begins.

The *Target HR* is set by the next parameter. The setting *Duration of training* controls the duration of the training phase (TP). The setting *Cool* controls the period of time in which the is reduced to 0 watts after termination of training phase. The following settings


for Po – no/yes, Ask for HR – no/yes and Ask for weight – no/yes can be viewed when a new training is started and edited. Editing these parameters can be required when different persons want to use the training program.

Training-Program			
Initial load: 50 Watt			
Duration A1: 1min			
Increase of load: 10 Watt/min			
Duration A2: 5 min			
Training HR: 130 bpm			
Duration training: 10 min			
Duration recovery: 2 min			
Ask for Po: aus			
Ask for HR: ein			
Ask for weight: aus			
cancel			modify

heart (TP),

down load the Ask

Starting the Training Program in the Menu **Start program** > **Training**.

ENGLISH

**Please Note:**

**Warming-up Phase 2:**

The training heart rate (target HR) should be reached during warming-up phase 2. However, A1 and A2 should together not exceed 5 to 8 minutes. This time is mainly determined by the *Po* value and the *increase in load*. From a physiological point of view, an *increase in load* by 10 watts per minute is recommended for a healthy person but it should not exceed 15 watts/minute even for well trained individuals. If the period of 5 to 8 minutes is not reached with the settings selected for the first training, the initial load (*Po*) should be adjusted.

**Training Phase:**

The effective training is in the training phase. It is important that the body, heart and circulation are trained but no overload is exerted. The pulse-steady-state program is a very easy means to achieve this. If the performance decreases by more than 15% within 20 minutes from the start of the training, the load is too high. In this case, reduce the target HR until an acceptable value is reached. If the performance decrease is less than 10%, the target HR can be increased.

**Adjustment:**

During the warming-up phase A1, the *Po* load can be adjusted in 5-watt steps using the + and - arrow keys.

During the training phase (TP), the target heart rate (target HR) can be adjusted using the arrow keys.

## **8 Maintenance and Fault Clearing**

### **8.1 Measurement Check**

The unit's measuring technology should be checked every 24 months. The following checks should be performed in the process:

1. Check of the ergometer's overall mechanical condition including linear motors and lifting spindles
2. Display check (contrast, lighting, etc.)
3. Check of correct ergometer rotational speed display
4. Check of the measurement sensor for brake performance
5. Check of mechanical power loss of the ergometer's drive system
6. Electrical safety check
7. Check of the blood pressure recorder's pressure measuring unit
8. Check for tightness of the pneumatic system
9. Check of the safety symbols and markings on the casing
10. Writing of an inspection report

These checks and any recalibrating work necessary should only be performed by authorised and trained personnel with the special tools required for this purpose.

### **8.2 Cleaning the Device**

The surface of the casing can be cleaned with a soft cloth that is dry or moist. Commercially available cleaning agents for household appliances can be used. The saddle should be cleaned with a leatherette cleaner.

It is imperative to make sure that no water penetrates the device.

The keyboard foil should never be cleaned with petrol, nitro cleaner or acetone.

### **8.3 Cleaning the Blood Pressure Cuff (250 SE)**

The blood pressure cuff consists of a waterproof plastic foil. It can be cleaned with soapy water and a cloth. The water temperature must not exceed 30 °C. The cuff should not be plunged into suds to clean, as the Velcro straps might mat. If it should nevertheless be required to wash the cuff in water, the microphone must first be removed and the air admission pipe must be closed.

## 8.4 Repairing Electromagnetic Disorders

The unit is only designed for operation in the following electromagnetic environment:

### Radio frequency emission according to CISPR 11, group 1, class B.

**Group 1** means that the ergometer exclusively uses HF energy for its internal function. This makes its HF emission very low and unlikely to disturb electronic devices in the vicinity.

**Class B** means that the ergometer is suitable for use in any facilities including residential areas, even if it is directly connected to the public mains that also supplies residential buildings.

The unit is resistant to jamming in an electromagnetic environment if the following prerequisites are met:

- The voltage corresponds to the typical business or hospital environment, and the humidity is at least 30%, especially if the floors are synthetic.

If any disorders should occur nevertheless, especially in the vicinity of devices labelled with the symbol  "non-ionic electromagnetic radiation", check the recommended minimal distance according to the following table. More information is given in the service manual.

<b>Recommended safety distances between portable and mobile HF telecommunication devices and the ergometer</b>			
The ergometer is designed for operation in an electromagnetic environment with controlled HF disturbance. The customer or user can help avoid electromagnetic disturbances by keeping the minimum distance between portable and mobile HF telecommunication devices (senders) and the ergometer, depending on the output performance of the communication device as indicated below.			
HF source	Rate [MHz]	Rated power P of the sender [W]	Distance [m]
Microcellular phone, CT1+, CT2, CT3, wireless headphones	885–887 MHz	0.01	0.23
Babyphone	27–41 MHz	0.1	0.37
Bluetooth systems (wireless mouse, wireless keyboard, handsfree set)	2400–2500	0.0025	0.38
Cordless DECT telephone, WLAN (laptop, PDAs), UMTS phone	1880–2500	0.25	1.17
Mobile phone	850/1900	1.2	1.8
Mobile phone, GSM850, NMT900, DCS 1800	850/900/1800	1	2.3
Mobile phone, GSM 900	900	2	3.3
Walkie-talkie (rescue services, police, fire brigade, maintenance services)	81–470	5	2.6
Mobile radio system (rescue services, police, fire brigade)	81–470	100	11.7

## 8.5 Checking and Setting the Supply Voltage

On delivery, the unit is set for the local supply voltage (110/115 ~ or 230/240 V ~). The current voltage setting is recorded on the mains module. The power unit covering lid on the bottom of the unit must be opened to convert the voltage. The voltage can then be changed on the power unit on-board with a special voltage selector switch.

Voltage conversion should only be performed by authorised and trained personnel.

## 8.6 Changing a Mains Fuse

The fuse switch is located in the centre of the mains module. The lid can be pried out of its lock-in position using a small screwdriver. It can then be pulled out of the fuse well. Two fuses are located in the fuse switch. After a continuity check, change the fuses if necessary. Return the fuse switch to the well and press it into the lock-in position.


Replace fuses only with other fuses of the same type with the same specification: 2 x 1.25 AT for 230 V, or 2 x 2.5 AT for 110 V.

## 8.7 Disposal Information

Devices that are no longer usable can be returned to ergosana for disposal. Alternatively, the device can be taken to an approved disposal location.

The device's control console contains a buffer battery, which must be disposed of separately.

## 9 Appendix

### 9.1 Technical Customer Service and Sales Locations

ergosana products are also sold as OEM products with other brand names. These devices are exclusively sold by authorised agents. The agents are trained in the service of our devices. Please contact one of these specialist dealers if your unit requires servicing.

If this is not possible please contact the company's central service department:

Service department  
ergosana GmbH  
Truchtelfinger Str. 17  
D-72475 Bitz

Phone: +49 74 31 9 89 75 13  
Fax: +49 74 31 9 89 75 15  
<http://www.ergosana.de/>

Notice