DEUTSCH

## Fahrradergometer Sana Bike 320 F

---

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ENGLISH

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**ENGLISH**

**CE Conformity Declaration**

ergosana GmbH herein declare that the medical products (Class IIa) of the ergometer system Sana Bike 320 F have been designed and manufactured in accordance with the relevant requirements of the EC directive 93/42/EEC Annex I.


This declaration loses its validity if the above devices are modified without ergosana's consent.

**The medical products are checked by the notified body DEKRA and they bear the CE mark CE 0124.**

Manufactured by:  
ergosana GmbH

Manufacturer's signature:

.....  
Harald Neukirchner  
Quality manager

Truchtelfinger Str. 17  
D-72475 Bitz

**Table of Contents**

**1 General.....4**

1.1 Intended use ..... 4

1.2 Contra-indications ..... 4

1.3 Features ..... 4

1.4 Instruction ..... 4

1.5 Maintenance ..... 4

**2 Product description .....5**

2.2 Accessories ..... 5

2.3 Potential Equalisation..... 5

2.4 Technical specifications ..... 6

2.5 Signs and symbols..... 6

**3 Installation.....7**

3.1 Location..... 7

3.2 Assembly instructions ..... 7

3.2.1 Unpacking and assembling..... 7

3.2.2 Adjusting the saddle and handle bars ..... 7

3.2.3 Connection..... 8

**4 Unit Components .....9**

4.1 Display ..... 9

4.2 Rotational speed display on control console ..... 9

4.2.1 Keys and Display..... 10

4.2.2 Setting the Language..... 10

**5 Safety Notes ..... 11**

5.1 Precautions during operation..... 11

5.2 Safety precautions when operating with other devices ..... 11

5.3 Maintenance safety precautions..... 11

5.4 Interference ..... 11

**6 Ergometry.....12**

6.1 Defining the Automatic Load Programs ..... 12

6.2 Recommended Settings..... 14

6.3 Remote Operation ..... 15

6.3.1 Explanation ..... 15

6.3.2 Prerequisites ..... 15

6.4 Training program (option) ..... 15

6.4.1 Training with constant heart rate (Pulse-Steady-State) ..... 15

6.4.2 Configuring a training program on the ergometer..... 16

**7 Maintenance and fault clearing.....19**

7.1 Measurement check..... 19

7.2 Cleaning the device ..... 19

7.3 Checking and setting the supply voltage ..... 19

7.4 Changing a mains fuse ..... 19

7.5 Eliminating electromagnetic interferences ..... 20

7.6 Disposal information ..... 21

**8 Appendix .....22**

8.1 Technical Customer Service and Sales Locations ..... 22

ENGLISH

## 1 General

The bicycle ergometer Sana Bike 320 F is a state of the art, high-performance ergometer. The device meets the highest quality standards in precise physical exertion tests.

### 1.1 Intended use

Ergometer Sana bikes are bicycle ergometers that are intended for defined exercise ergometry during a patient's examination and therapy. These products are used in practices, clinics, therapy and rehabilitation centres. The ergometers are operated by physicians and medical personnel.

### 1.2 Contra-indications

**In the event of the following contra-indications, NO exercise test must be performed:** in the case of acute cardiac infarction or unstable angina pectoris, serious hypertonia at rest, carditis, insufficiency of the heart, serious valvular heart defect, serious cardiac arrhythmia at rest, aortic aneurysm or other manifest cardiovascular diseases.

### 1.3 Features

**The following characteristics make the unit exceptional:**

- Attractive design
- Comfortable mounting and dismounting
- Stable steel construction, compact drive unit
- Stable clamps on saddle and handle bars
- Standardised saddle pipe (change of saddle possible at any time)
- Impact- and scratch-proof casing of glass fibre reinforced plastic with acrylic surface, easy to clean
- New, high-performance control electronics
- Graphic display featuring visual representation of ergometry data
- Easy operation via menu guidance
- Remote operation – personalised programs – training programs
- Performance range from 1 to 400
- Guaranteed accuracy (error < 3%, not less than 3 watts in the independent rpm range)
- Very quiet drive mechanism
- Pleasant pedalling sensation due to large gyrating mass
- Galvanically isolated RS-232 interface, secure data transfer

### 1.4 Instruction

Before the initial operation, carefully read through this user guide, paying special attention to the warnings and safety instructions.

### 1.5 Maintenance

This is a low-maintenance device. You will find detailed maintenance instructions in section 7.

## 2 Product description

### 2.1 Unit components

1. Handle bar
2. Saddle
3. Clamp for saddle height adjustment
4. Mains connector, potential equalization, RS-232
5. Base adjustor
6. Running gear with lockable heavy rollers
7. Handle bar clamp
8. Control panel with LCD display and foil keyboard with control elements and LED for rotational speed display

### 2.2 Accessories

#### Every device comes with:

- Power cable with European plug
- User guide
- Inspection report

### 2.3 Potential Equalisation

A standard potential equalisation stud is located on the rear panel, next to the power connection unit. It is marked with a green/yellow information sign. Using an earthing cable, the ergometer can be connected to the potential equalisation of the examining room, which serves as a common earthing point for all other mains-operated devices in the room to ensure that all devices have the same earthing potential.

## 2.4 Technical specifications

### Bicycle ergometer with blood pressure measurement according to DIN VDE 0750-238

<b>Braking principle</b>	Computer-controlled brakes with permanent measurement of torque. Braking performance is independent of revolutions per minute. (see load range precision)
<b>Power range</b>	1 till 400 Watts
<b>Load range</b>	independent rpm range 20 to 400 Watts
<b>Range of revolutions</b>	30 to 130 /min
<b>Load precision</b>	3%, not less than 3 watts (in the independent rpm range)
<b>Load parameters</b>	<ol style="list-style-type: none"> <li>1. In keeping with set internal load program</li> <li>2. Parameters from external master unit over interface, in 1 watt steps.</li> <li>3. Manual in 5-watt and 25-watt steps</li> </ol>
<b>Load software</b>	5 freely programmable ergometry programs
<b>Time intervals</b>	1 min to 99 min
<b>Display</b>	LCD display with graphic capability, 320 x 240 pixels, CCFT background lighting.
<b>Pulse measurement</b>	Option with pulse belt, display accuracy +/- 5%
<b>Adjustable seat and handle bar</b>	Infinitely adjustable for heights ranging from 150 cm to 210 cm
<b>Long-term accuracy</b>	Continuous torque control and equalisation according to weight
<b>Power supply</b>	230 VAC 50-60 Hz , 115 VAC 50-60 Hz The unit is suitable for use in networks according to CISPR 11, group 1, class B.
<b>Electric inputs/outputs</b>	RS-232 (galvanically isolated), analogue remote boot
<b>Base dimensions</b>	40 x 83 cm
<b>Weight</b>	43 kg

## 2.5 Signs and symbols

In this section, the signs and symbols used in connection with this device are explained:

Mains operated, alternating current

Potential equalisation connection (earth)

BF classified component

Warning! Follow the instructions in the documentation.

93/42/EEC for medical products 0124  
CE 0124 DEKRA

**IPX0**

Protection class of the casing: IPX0

## 3 Installation

### 3.1 Location

Install the device in a suitable position (refer to safety instructions in section 5).

The unit must not be stored or operated in wet, moist or dusty surroundings, nor should the unit be exposed to direct sunlight or other sources of warmth.

The unit must not come into contact with acidic vapours or fluids.

The unit should not be placed near X-ray units, large transformers or electrical motors. There must be a distance of at least one meter between the unit and the mains network.

### 3.2 Assembly instructions

#### 3.2.1 Unpacking and assembling

Install the control panel after unpacking the unit. In order to do so, fix the back of the control panel over the handle bar pipe, insert the two tabs on the back of the control panel into the handle bar pipe and press them downward to their stop. Then connect the connection cables with the control panel and screw on the cover.

Note:

Mount the control panel with the operator's side facing the front so the keys and display are easily accessible by the person operating the machine. If the device is used as a training ergometer, the panel can be mounted with the operator's side facing the rider.

Screw the handle bar on with a rail clip and a hex socket screw. Make sure that the slit created by screwing the pieces together is the same width on the bottom as at the top due to the clamping action. Then screw on the clamp and fix the handle bar in place. Adjust the clamp lever so that it faces downward.

Pull the saddle to normal height, adjust the clamp lever downwards after fixing the saddle in place.

With the help of the base adjustor on the rear lower side of the ergometer, adjust the device so that there is no gap between it and the floor. The ergometer is then fully stabilised.

#### 3.2.2 Adjusting the saddle and handle bars

The saddle can be adjusted upwards or downwards for optimal exercise position. The handle bar can be rotated by 360 degrees. This allows the optimal adjustment for exercise for persons ranging from 150 cm to 210 cm.

The handle bar is adjusted by means of the clamp situated directly below the handle bar. The saddle is adjusted using a clamp protruding from the casing below the saddle rod. After loosening the respective clamp, the handle bar can be adjusted to a suitable position and the saddle pipe can be moved up or down to the desired height. After adjustment, tighten the clamps again.

The saddle clamp is designed for patients weighing up to 160 kilograms. To attain optimal clamping action, only moderate strength is needed to tighten the clamps.

It is recommended to position the clamp in clamped position with the lever positioned vertically downward. This is achieved by pulling the clamp out of its lock-in position and then turning it until the lever is positioned vertically downward. Secure clamping is ensured when the lever is once again turned to this position during the next adjustment procedure.

### 3.2.3 Connection

Establish potential equalization (refer to section 2.3) and plug the supplied power cable into an earthed socket. As the device is preset to the local mains voltage (refer to section 7.3), you can switch it on using the main switch on the rear side. The ergometer is now ready for operation.

## 4 Unit Components

### 4.1 Display

The console is facing the user (not the patient) (factory setting)

It is possible to turn the console by 180 degrees for special applications such as patient training, etc., so that the patient can reach the control elements and see the display.

All of the control electronics for all operation modes of the ergometer are located in the control console.


A backlit graphic LCD display with 320 x 240 pixels and a surface of 100 x 75 mm featuring the complete range of available information is located on the front side under a foil keyboard with a see-through window. The control elements for the adjustment and operation of the ergometer are located on the foil keyboard.

A clear LCD display with large figures is located on the control console's upper side. It shows the rider the current number of pedal rotations per minute.

### 4.2 Rotational speed display on control console

$n$  = crank rotations per min.

$U \text{ min}^{-1}$

### 4.2.1 Keys and Display


The display shows all operating and measuring procedures. When the unit is switched on, the basic menu is displayed. The desired menu option is selected using the arrow keys.

The right and left buttons on the control panel have been arranged so that the function executed by the corresponding key is shown directly above the key on the lower side of the display. The desired menu option is opened by pressing the respective key. All menu options required for the correct operation of the device are located in the "Settings" menu. The settings in this menu should only be defined by service technicians and trained staff. All further setting options are described in section 6.

During an exercise recording, the current load curve is displayed (trend graph). If the heart rate is loaded (optional), it will also be displayed with a different line width than the load curve. This makes the monitoring of the load parameters during a stress test particularly easy.

The upper part of the display is alphanumeric. This part contains several boxes indicating the current speed, the current load, the running period and, if the measurement option is installed, the current heart rate.

It is possible to change the display mode from Trend graph to Numerical. See page 18 Menu Settings/Display mode.

- 1 = Confirmation key
- 2 = "Up" and "down" cursor keys
- 3 = LCD display

### 4.2.2 Setting the Language

The preset menu language is "deutsch".

To change the language, highlight the "Einstellungen" menu option using the "ab" arrow, and confirm with "OK".

In the "Einstellungen" menu, select the "Sprache" option and confirm with "OK". Then select the desired language using the "auf" or "ab" arrow, and confirm with "OK".

The menu is now displayed in the selected language.

## 5 Safety Notes

### 5.1 Precautions during operation

Before using the unit, make sure that the Medical Product Consultant has conducted an introduction in regard to function and safety precautions.

The unit should not be used if there are any doubts as to its being isolated from earth or the suitability of the power cable.

The supplied power cable meets valid regulations for medical applications.

The unit is not intended for use in wet rooms, outdoors or in areas where there is danger of explosion.

Before the initial operation, the device must be adjusted using the base adjusters at the back to grant absolute stability.

When the saddle is exchanged, make sure that the screws at the saddle are tightened hard enough so that the saddle cannot be moved on the saddle pipe.

To move the handle bar and the saddle, loosen the clamps and retighten them afterwards. It is recommended to position the clamps with the levers positioned vertically downward. Secure clamping is ensured when the handles are then again turned to this position during the adjustment procedure after every movement.

The holding straps on the pedals must fit perfectly across the upper side of the shoe and be fastened with a Velcro strap.

### 5.2 Safety precautions when operating with other devices

When several devices are coupled, there is a risk that the leakage currents may add up.

The RS-232 interface, which can be used for communication with other devices, is isolated to ensure the patient's safety.

External devices may only be connected with the interface cables supplied by ergosana.

Portable communication devices, HF radios and devices labelled with the symbol  (non-ionic electromagnetic radiation) can affect the operation of this device (see section 7.5).

### 5.3 Maintenance safety precautions

The device must be turned off and the power plug disconnected before cleaning with liquid cleaning agents.

Only use standard cleaning agents for plastic surfaces.

The unit may only be opened, repaired and serviced by authorised and trained personnel.

### 5.4 Interference


The unit meets EMC regulations for medical products to ensure protection against emission and radiation. Special caution should be taken when using this unit in combination with high-frequency devices.

## 6 Ergometry

This section describes ergometry with the internal ergometer software.

### 6.1 Defining the Automatic Load Programs


The following text field appears in the ergometer after it is switched on:

The menu item "Start ergometry" is highlighted with a black bar. This means the item is activated. It is possible to select an exercise program directly from here and begin ergometry.


Beforehand, however, the five varying exercise programs must be defined according to the examiner's needs and wishes. On delivery of the unit, each program is set with normal values which cannot harm the patient in case they are accidentally activated.

Activate "Define program" via the arrow keys  .


Confirm with the right "ok" button.

The program selection menu appears.

The exercise programs 1–5 are set according to the parameters of standardised ergometry or according to individual parameters in such a manner that the correct programs for the different patient groups can be called up at the touch of a button.

Confirm program 1 with the right button “ok”, and the window containing the settings menu for ergometry program no. 1 appears.

Check the existing settings for:

- Type of load (Stage or Ramp)
- Initial load
- Increase in load
- Stage time
- Recovery load

Push the “change” button if the settings need to be altered.

The individual parameters then appear one after the other. Select the desired values with the arrow buttons and confirm by pressing the “ok” button. The next parameter then appears. The settings menu can be exited at any time by pressing the “cancel” button.

## 6.2 Recommended Settings

An unlimited number of values can be saved in the programs 1 to 5, as shown in the examples below:

Program number	Initial load [watts]	Increase in load [watts]	Stage time [min]	Recovery load [watts]
1	30	10	1	20
2	25	25	2	25
3	50	25	2	25
4	50	50	3	50
5	75	50	3	50

The program always returns to the initial configuration after the individual programs have been set.

“Start ergometry” is highlighted black on the LCD display. An ergometry exercise can be started immediately by pressing the “ok” button.

## 6.3 Remote Operation

### 6.3.1 Explanation

Remote operation means that the ergometer is externally controlled via the digital RS-232 interface or USB connector, i.e. all commands for the level of exertion are transmitted from a separate "master unit".

This type of operation is used primarily when several ergometers are controlled from a central workstation, the ECG device has its own ergometry-control program and the ergometer and ECG unit – and possibly other devices such as spirometry devices – can be combined to form an ergometry or spirometry station.

### 6.3.2 Prerequisites

When the remote operation mode is selected, the unit used (an ECG device or a PC) must be connected with the ergometer via an interface cable. For our ergometers, this is a type RS-232 interface, which is isolated to ensure the patient's safety.

The appropriate baud rate must be selected in the program "Setting" under "Interface". The appropriate transmission protocol must then be selected in the same menu under "Instruction Set". The P 10 setting contains the operation mode for SCHILLER devices, among others. Our own transmission protocol is located under "ergosana".

Information on interface and instruction set can be obtained from the operating data of the master unit.

When the settings have been correctly defined, the ergometer automatically switches to "remote operation" when the first command comes through over the interface. The ergometry screen is opened in the display, showing current load and possibly pulse rate data in both alphanumeric and graphic forms. The execution of the internal programs is disabled in this operating mode.

Remote operating mode is ended by either pressing the command button "End" or by switching the device off.

## 6.4 Training program (option)

The use of the training program requires a pulse signal receiver (Polar system), which is integrated in the control panel. The receiver can be ordered together with the device or added later on.

The patient wears a transmitter belt, which is positioned on the skin below the chest. The effective radius of the signals transmitted by the belt and received in the control panel is approx. 70 cm. Please note that if the skin is dry, contact problems may occur between the belt and the skin at the beginning of the training. Therefore, if the pulse transmission is unsteady or disturbed, moisten the contact surfaces of the transmitter belt with contact spray or water.

### 6.4.1 Training with constant heart rate (Pulse-Steady-State)

A bicycle ergometer training with constant heart rate (pulse-steady-state method) in the individual training range is a highly efficient and risk-free training method for the cardiopulmonary system. Ask your physician which is the suitable training heart rate for you. He or she will determine it for you using an exercise test.


## 6.5 Configuring a training program on the ergometer

To define a training program, select the menu option **Define program**. Then select the menu option **Training** using the arrow keys. The below settings panel is displayed. To set or change the parameters in order, press **modify**.

The first setting defines the *Initial load*. e.g. 50 watts. The setting *Duration A1* controls the duration of warming-up phase 1, e.g. 1 min. The setting *Increase in load* controls by how many watts per minute the load is increased. This is warming-up phase 2. Using the next setting, *Duration A2*, this phase can be limited. By default, the warming-up phase 2 is finished when the target HR is reached. When the target HR is reached, the load applied is automatically reduced by 10%. This prevents an "excessive increase in heart rate" during the training phase (TP) which now begins. The *Target HR* is set by the next parameter. The setting *Duration of training* controls the duration of the training phase (TP). The setting *Cool down* controls the period of time in which the load is reduced to 0 watts after termination of the training phase. The following settings Ask for Po – no/yes, Ask for HR – no/yes and Ask for weight – no/yes can be viewed when a new training is started and edited. Editing these parameters can be required when different persons want to use the training program.

Training-Program			
Initial load: 50 Watt			
Duration A1: 1min			
Increase of load: 10 Watt/min			
Duration A2: 5 min			
Training HR: 130 bpm			
Duration training: 10 min			
Duration recovery: 2 min			
Ask for Po: aus			
Ask for HR: ein			
Ask for weight: aus			
<b>cancel</b>			<b>modify</b>

Starting the Training Program in the menu **Start program** > **Training**.

ENGLISH

**Please note:**

**Warming-up phase 2:**

The training heart rate (target HR) should be reached during warming-up phase 2. However, A1 and A2 should together not exceed 5 to 8 minutes. This time is mainly determined by the *Po* value and the *increase in load*. From a physiological point of view, an *increase in load* by 10 watts per minute is recommended for a healthy person but it should not exceed 15 watts/minute even for well trained individuals. If the period of 5 to 8 minutes is not reached with the settings selected for the first training, the initial load (*Po*) should be adjusted.

**Training phase:**


The effective training is in the training phase. It is important that the body, heart and circulation are trained but no overload is exerted. The pulse-steady-state program is a very easy means to achieve this. If the performance decreases by more than 15% within 20 minutes from the start of the training, the load is too high. In this case, reduce the target HR until an acceptable value is reached. If the performance decrease is less than 10%, the target HR can be increased.

**Adjustment:**

During the warming-up phase A1, the *Po* load can be adjusted in 5-watt steps using the + and - arrow keys.

During the training phase (TP), the target heart rate (target HR) can be adjusted using the arrow keys.

Select the menu "Settings" and then "Display mode" to select either trend graphs or numerical presentation. With the numerical presentation, individual values are easier to discern.

**Numerical display mode**

Program: 1

Duration <b>10</b>	Speed <b>56</b>
Load [watts] <b>80</b>	Pulse [1/min] <b>95</b>
Energy [kcal] <b>86</b>	
Systole [mmHg] <b>145</b>	Diastole [mmHg] <b>84</b>
Cancel	Start

## 7 Maintenance and fault clearing

### 7.1 Measurement check

The unit's measuring technology should be checked every 24 months. The following checks should be performed in the process:

1. Check of the overall mechanical condition of the ergometer
2. Display check (contrast, lighting, etc.)
3. Check of correct ergometer rotational speed display
4. Check of the measurement sensor for brake performance
5. Check of mechanical power loss of the ergometer's drive system
6. Electrical safety check
7. Check of the safety symbols and markings on the casing
8. Writing of an inspection report

+ These checks and any recalibrating work necessary should only be performed by authorised and trained personnel with the special tools required for this purpose.

### 7.2 Cleaning the device

The surface of the casing can be cleaned with a soft cloth that is dry or moist. Commercially available cleaning agents for household appliances can be used. The saddle and the handlebar should be cleaned with a leatherette cleaner.

+ It is imperative to make sure that no water penetrates the device.

+ The device should never be cleaned with petrol, nitro cleaner or acetone.

### 7.3 Checking and setting the supply voltage

On delivery, the unit is set for the local supply voltage (110/115 ~ or 230/240 V ~). The current voltage setting is recorded on the mains module. The power unit covering lid on the bottom of the unit must be opened to convert the voltage. The voltage can then be changed on the power unit on-board with a special voltage selector switch.

+ Voltage conversion should only be performed by authorised and trained personnel.

### 7.4 Changing a mains fuse

The fuse switch is located in the centre of the mains module. The lid can be prised out of its lock-in position with the help of a small screwdriver. It can then be pulled out of the fuse well. Two fuses are located in the fuse switch. After a continuity check, change the fuses if necessary. Return the fuse switch to the well and press it into the lock-in position. If the fuse won't stay in position, immediately inform the after-sales service.

+ Replace fuses only with other fuses of the same type with the same specification:  
2 x 1.25 AT for 230 V, or 2 x 2.5 AT for 110 V.

## 7.5 Eliminating electromagnetic interferences

The unit is only designed for operation in the following electromagnetic environment:

### Radio frequency emission according to CISPR 11, group 1, class B.

**Group 1** means that the ergometer exclusively uses HF energy for its internal function. This makes its HF emission very low and unlikely to disturb electronic devices in the vicinity.

**Class B** means that the ergometer is suitable for use in any facilities including residential areas, even if it is directly connected to the public mains that also supplies residential buildings.

The unit is resistant to jamming in an electromagnetic environment if the following prerequisites are met:

- The voltage corresponds to the typical business or hospital environment, and the humidity is at least 30%, especially if the floors are synthetic.

If any disorders should occur nevertheless, especially in the vicinity of devices labelled with the symbol  "non-ionic electromagnetic radiation", check the recommended minimal distance according to the following table. More information is given in the service manual.

<b>Recommended safety distances between portable and mobile HF telecommunication devices and the Sana Bike</b>			
The Sana Bike ergometer is designed for operation in an electromagnetic environment with controlled HF disturbance. The customer or user can help avoid electromagnetic disturbances by keeping the minimum distance between portable and mobile HF telecommunication devices (senders) and the ergometer, depending on the output performance of the communication device as indicated below.			
<b>HF source</b>	<b>Rate [MHz]</b>	<b>Rated power P of the sender [W]</b>	<b>Distance [m]</b>
Microcellular phone, CT1+, CT2, CT3, wireless headphones	885–887 MHz	0,01	0,23
Babyphone	27–41 MHz	0.1	0.37
Bluetooth systems (wireless mouse, wireless keyboard, handsfree set)	2400-2500	0.0025	0.38
Cordless DECT telephone, WLAN (laptop, PDAs), UMTS phone	1880-2500	0.25	1.17
Mobile phone USA	850/1900	1,2	1.8
Mobile phone, GSM850, NMT900, DCS 1800	850/900/1800	1	2.3
Mobile phone, GSM 900	900	2	3.3
Walkie-talkie (rescue services, police, fire brigade, maintenance services)	81-470	5	2.6
Mobile telephone system (rescue service, police, fire brigade)	81-470	100	11.7

---


## 7.6 Disposal information

Devices that are no longer usable can be returned to ergosana for disposal. Alternatively, the device can be taken to an approved disposal location.

- + The control console of the device contains a buffer battery, which must be disposed of separately.

## 8 Appendix

The family of characteristics for the area of the braking moment control for the Sana Bike ergometer type 400 Watt

### 8.1 Technical Customer Service and Sales Locations

ergosana products are also sold as OEM products with other brand names. These devices are sold by authorised agents. The agents are trained in the service of our devices. Please contact one of these specialist dealers if your unit requires servicing.

If this is not possible, please contact the company's central service department:

Service department  
ergosana GmbH  
Truchtelfinger Str. 17  
D-72475 Bitz

Phone +49 74 31 9 89 75 0  
Fax +49 74 31 9 89 75 15  
<http://www.ergosana.de/>  
[info@ergosana.de](mailto:info@ergosana.de)